


THE HAMMOCK







The Hammock is the quick and easy way to carry a Tourniquet or other similar sized objects, on an existing belt or MOLLE mounted pouch. The Hammock requires no more additional mounting space and protects the tourniquet from the environment. With the Tourniquet Hammock, you can keep your tourniquet readily accessible right below your Micro Trauma Kit or other Molle mounted equipment. As an added benefit, the Tourniquet Hammock keeps the tourniquet better protected against damaging UV rays, weather, and dirt and debris.

Weight	19g
Part #	HW-TQ-01-XX
Size (empty)	5" x 3" x 3/8"
Size (full)	5" x 3" x 2.5"

Ready to attach to any Molle or Belt attachment

Designed for Rapid Deployment

Fits CAT tourniquet & other similar sized objects

Solution Dyed CORDURA®

ULTRAcomp®



Medical







